

Perma Evolution

Description

This permanent root filling material combines epoxy-amine technology with innovative benefits for the dental professional.

Application

After preparing the root canal using standard methods, and if necessary using a temporary filling material (e.g. Asphalene temp), the root canal is filled with Perma Evolution and guttapercha points. Before use, the cap of the double-chamber syringe must be removed by turning counter-clockwise. A mixing cannula supplied with the product is attached by inserting the spring of the mixing cannula into the groove in the syringe. The cannula is tightened by turning in a clockwise direction. When you press the piston, both pastes are mixed in the mixing cannula in a ratio of 1:1 to form a ready-to-use paste. The paste can be applied and worked from a mixing block using a lentulo or guttapercha point. Alternatively, the material can be applied directly to the canal openings and coronal part of the canal using Endo tips. After treatment is completed, the cap must be screwed back on to the syringe the mixing cannula discarded.

Properties

- free of iodoform and paraformaldehyde
- excellent biocompatibility
- microbicidal action
- outstanding dimensional stability
- strong adhesion to the dentin
- radiopaque
- easy retreatment
- suitable for all orthograde root filling techniques
- fulfills the requirements of EN ISO 6876 for dental root canal sealing materials
- working time: 15 minutes
- setting time: 24 hours



Composition

Paste 1 :

4-[2-(4-hydroxyphenyl)propan-2-yl]phenol-Epichlorohydrinresin, Alkylglycidyl ether, Barium sulfate, tri-calcium phosphate, diphenylolpropan-diglycidyl ether

Paste 2:

Polyalkoxyalkylamine-copolymer, 5-amino-1,3,3-trimethylcyclohexanmethylamin, Aqua, Barium sulfate, tri-calciumphosphate, nanodispers silicon dioxide, polyhexamethylene biguanides-hydrochloride

Form of Delivery

REF 555 Automix syringe with 5 ml content, 20 mixing cannulas and 20 endo tips

REF 555.0 Refill with automix syringe with 5 ml content, without cannulas and endo tips

Accessories

REF 555.1 Refill with 25 mixing cannulas

REF 555.2 Refill with 25 endo tips

Shelf life

24 months

Classification

Medical device, class IIa, Regel 8,
CE Certificate: G1 17 09 49628 010,
TÜV SÜD Product Service GmbH, München

Manufacturer

Alfred Becht GmbH
Carl-Zeiss-Str. 16
D-77656 Offenburg
Tel.: +49 781 60586-0
Fax.: +49 781 60586-40
Germany
www.becht-online.de

Starter-Kit

Automix syringe with 3g content, 2 mixing cannulas and 2 endo tips